
Midmark M11 - v1.0.0
Total Cycles: 597

Name: MIDMARK
Sterilizer ID: V2178198
01/12/2021 10:22 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 136.8 F 0.0
2:00 139.0 F 1.2
4:00 159.2 F 3.8
6:00 186.1 F 4.6
8:00 209.4 F 3.2
10:00 219.7 F 5.9
12:00 228.0 F 6.0
14:00 227.1 F 5.0
16:00 241.2 F 10.8
18:00 252.7 F 16.3

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 15.9
1:00 252.2 F 16.0
2:00 252.7 F 16.4
3:00 251.6 F 15.8
4:00 252.7 F 16.3
5:00 251.6 F 15.8
6:00 252.7 F 16.4
7:00 251.8 F 15.9
8:00 252.6 F 16.5
9:00 252.3 F 16.1
10:00 251.6 F 15.9
11:00 252.7 F 16.4
12:00 252.2 F 16.1
13:00 251.6 F 15.8
14:00 252.7 F 16.5
15:00 252.5 F 16.3
16:00 251.9 F 15.9
17:00 252.0 F 16.3
18:00 252.8 F 16.4
19:00 252.3 F 16.1
20:00 251.7 F 15.8
21:00 252.4 F 16.5
22:00 252.7 F 16.4
23:00 252.2 F 16.1
24:00 251.6 F 15.8
25:00 252.6 F 16.5
26:00 252.7 F 16.4
27:00 252.2 F 16.1
28:00 251.7 F 15.8
29:00 252.5 F 16.5
30:00 252.7 F 16.4

Min 251.5 F 15.8
Max 252.9 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 172.8 F 0.0
5:00 188.7 F 0.0
10:00 180.2 F 0.0
15:00 202.1 F 0.0
20:00 203.6 F 0.0
25:00 208.4 F 0.0
30:00 206.3 F 0.0
Drying Complete

FILLING: 0:53

HEATING: 18:17
STERILIZING: 30:00
VENTING: 2:50
DRYING: 30:00
TOTAL CYCLE: 01:22:00
CYCLE COMPLETE

01/12/2021 11:44 AM
ACCEPTED: X
REJECTED:
Operator: U01