
Midmark M11 - v1.0.5
Total Cycles: 3324

Name: RITTER2
Sterilizer ID: V2313228
01/26/2024 09:34 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 74.5 F 0.0
2:00 88.1 F 0.0
4:00 123.8 F 0.2
6:00 118.2 F 1.7
8:00 145.8 F 4.6
10:00 175.4 F 6.7
12:00 209.1 F 2.5
14:00 220.2 F 6.0
16:00 217.1 F 2.6
18:00 228.9 F 6.7
20:00 226.6 F 5.1
22:00 238.6 F 10.0
24:00 249.3 F 15.1

STERILIZING
mm:ss Degrees PSI
0:00 251.6 F 16.6
1:00 252.0 F 16.5
2:00 252.6 F 16.8
3:00 252.6 F 17.0
4:00 252.0 F 17.1
5:00 251.5 F 16.5
6:00 252.7 F 16.9
7:00 252.2 F 17.1
8:00 251.9 F 16.6
9:00 252.3 F 17.0
10:00 251.5 F 16.6
11:00 252.6 F 17.0
12:00 251.4 F 16.6
13:00 252.2 F 17.0
14:00 251.4 F 16.6
15:00 252.2 F 17.0
16:00 251.7 F 16.6
17:00 252.3 F 17.1
18:00 251.8 F 16.7
19:00 252.0 F 17.1
20:00 252.5 F 16.8
21:00 251.5 F 16.8
22:00 253.0 F 17.0
23:00 251.7 F 16.6
24:00 252.3 F 17.1
25:00 252.5 F 16.8
26:00 251.4 F 16.7
27:00 252.6 F 17.0
28:00 252.1 F 16.7
29:00 251.9 F 17.1
30:00 252.7 F 16.9

Min 251.2 F 16.5
Max 253.0 F 17.1

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 187.5 F 0.0
5:00 214.8 F 0.0
10:00 218.8 F 0.0
15:00 197.5 F 0.0
20:00 195.0 F 0.0
25:00 193.2 F 0.0
30:00 208.4 F 0.0

Drying Complete

FILLING: 0:49
HEATING: 25:40
STERILIZING: 30:00
VENTING: 2:59
DRYING: 30:00
TOTAL CYCLE: 01:29:28
CYCLE COMPLETE

01/26/2024 11:03 AM
ACCEPTED:
REJECTED:
Reason: _____
Operator: _____